

# CSM/CVSM SERIES DEWARs

## SUPERCONDUCTING MAGNET DEWARs


Go To:  
[CVSM](#)

### CSM


#### CSM LN<sub>2</sub> SHIELDED DEWARs




Model	"A" Inside Diameter (Inches)	"B" Outside Diameter (Inches)	"C" Inner Depth (Inches)	"D" Overall Height (Inches)
CSM648	6.00	12.00	48.00	54.00
CSM7548	7.50	13.00	48.00	54.00
CSM8548	8.50	14.00	48.00	54.00
CSM9548	9.50	16.00	48.00	54.00
CSM10548	10.50	16.00	48.00	54.00
CSM12548	12.50	18.00	48.00	54.00
CSM1548	15.00	20.00	48.00	54.00



MAGNET DEWARs



MAGNET DEWAR, LN<sub>2</sub> SHIELDED WITH BELLY



MAGNET DEWAR, LN<sub>2</sub> SHIELDED WITH BELLY AND TAIL

[MORE OPTIONS AVAILABLE](#)