

# CLPB 15 LN<sub>2</sub> Container




- Ideal for conventional straight liquid dispensing
- Sizes range from 15 L to 500 L
- Convenient fill/withdrawal valve
- Vent valve with full tricock prevents overfilling
- Pressure gauge and relief devices for safety
- CLPB series adds a self-pressurizing system for convenience and transfer efficiency
- CLPB series feature adjustable regulator to maintain 10–25 PSI
- Dewar can be customized to fit your specific needs

[View Another Size?](#)

- [25 L](#)
- [50 L](#)
- [80 L](#)
- [110 L](#)
- [160 L](#)
- [240 L](#)
- [500 L](#)

	CLPB
Net Liter Capacity	15 L
Gross Liter Capacity	17 L
Outside Diameter	12"
Height	30"
Empty Weight	37 lbs
Weight Full LN <sub>2</sub>	63.7 lbs
Weight Full LOX	75 lbs
Weight Full LARG	82 lbs
N.E.R. (Normal Evaporation Rate)	1.4 L/day
M.A.W.P. (Maximum Allowable Working Pressure)	25 PSI
Regulator Setting	22 PSI
Secondary Relief	50 PSI
Level Gauge	No
Liquid Fitting	—CGA 295 1/2" 45° Male Flare—
Vent Fitting	—3/8"NPT—
Material of Construction	—304 Grade Stainless Steel—

## Available Options:



Symbol	Description
V-1	Liquid Fill/Withdrawal
V-2	Fill/Withdrawal Valve
V-3	Pressure Building
V-4	Vent Valve
RD-1	Rupture Disk
RV-1	Relief Device
Pg-1	Pressure Gauge
PCV-1	Pressure Building Regulator
PBC-1	Pressure Building Coil

- Gravity Feed Port
- Top Vertical Withdrawal
- Electric Pressure Builder