

CH 600 LN₂, LOX, LARG Container

- Horizontal configuration allows large volumes of liquid to be portable and accessible
- Lower profile with 2" bottom clearance
- Conveniently located 1.5" o.d. neck opening with flanged termination
- Casters for easy mobility
- Unique support system insures durability without sacrificing thermal efficiency
- Lowest possible head leak and loss rates
- Stainless steel inner and outer jacket
- Maximum working pressure 25 PSI
- Sizes range from 300 to 2000 liters, and can be customized to suit your specific needs

**View
Another
Size?**

[300L](#)

[400L](#)

[500L](#)

[800L](#)

[1000L](#)

[1500L](#)

[2000L](#)

Net Capacity	600 L
Gross Capacity	660 L
Outside Diameter	32"
Overall Height	39"
Overall Length	82"
Empty Weight	697 lbs
Weight Full LN ₂	1663 lbs
Weight Full LOX	2105 lbs
Weight Full LARG	2395 lbs
N.E.R.	1.25%
M.A.W.P.	25 PSI
Safety Relief	100 PSI
Pressure Builder	20 PSI
Castor Size	5"
Liquid Fitting	CGA 295 1/2" 45 degrees Male Flare
Vent Fitting	3/8" NPT

Available Options:



- Removable Manifold
- Self-Pressurizing System
- Front-Mounted Plumbing
- Skid Mounting
- Differential Pressure Liquid Level Gauge
- Extra Access Port
- Top-Mounted Plumbing