

# CSM/CVSM SERIES DEWARS

## SUPERCONDUCTING MAGNET DEWARS

More Sizes Available


### CSM

U171006


#### CSM LN<sub>2</sub> SHIELDED DEWARS




Model	Inside Diameter (Inches)	Outside Diameter (Inches)	Inner Depth (Inches)	Overall Height (Inches)
CSM648	6.00	12.00	48.00	54.00
CSM7548	7.50	13.00	48.00	54.00
CSM8548	8.50	14.00	48.00	54.00
CSM9548	9.50	16.00	48.00	54.00
CSM10548	10.50	16.00	48.00	54.00
CSM12548	12.50	18.00	48.00	54.00
CSM1548	15.00	20.00	48.00	54.00



MAGNET DEWARS



MAGNET DEWAR, LN<sub>2</sub> SHIELDED WITH BELLY



MAGNET DEWAR, LN<sub>2</sub> SHIELDED WITH BELLY AND TAIL